



ETEST® Ceftazidime/Avibactam (0.016-256 µg/mL)

FDA 510(k) CLEARED

ETEST® CZA (0.016-256 µg/mL)



REF. 419046

Single pack | MIC range/concentration

30 Strips | Ceftazidime: 0.016-256 µg/mL
Avibactam: 4 µg/mL

Available from bioMérieux, Inc.

For determining ceftazidime/avibactam MIC for Enterobacteriaceae and *Pseudomonas aeruginosa*.



For more information, contact your local representative.

Customer Service: (800) 682-2666 | www.biomerieux-usa.com/etestcza



bioMérieux, Inc.
100 Rodolphe Street
Durham, NC 27712

ETEST® Ceftazidime / Avibactam

FDA 510(k) cleared!

REF# 419046 ETEST CZA 0.016-256 µg/mL.

Performance characteristics are provided below.

Organism	Strains	% Essential Agreement	% Category Agreement
	(N)	(EA)	(CA)
Enterobacteriaceae	1087	99.1	99.6
<i>Pseudomonas aeruginosa</i>	152	99.3	99.3

Note: EA = % of MIC values \pm 1 dilution of the reference method.

For more information, contact your local representative.

Customer Service: (800) 682-2666
www.biomerieux-usa.com/etestcza

bioMérieux, Inc.
100 Rodolphe Street • Durham, NC 27712 USA
Tel: (800) 682 2666 • Fax: (800) 968 9494
www.biomerieux-usa.com